3M Paradigm™ Temporization Material

Temporization Material

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ENGLISH

Product Description

Paradigm™ Temporization Material is a composite material for the corresponding temporary restorations. Based on multi-functional methacrylate esters, this two-component system has filling material-like properties and is available in the following colors: (e.g., A2, White) of the Garant™ Dispenser, the material can be applied directly from the cartridge without causing any bubbles to form. Paradigm™ Temporization Material temporary restorations are individually used as vacuum-formed or vacuum-cured composites (e.g., from the Filtek™ Flowable Composite product line).

For details on all mentioned products, please refer to the corresponding information for Use, Instructions for Use or the corresponding technical data sheets. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. The use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If the product contacts the eyes, remove and discard the gloves, wash the hands immediately with soap and water, and then re-glove. If allergic reaction occurs, seek medical attention as needed.

3M MSDS can be obtained at www.3m.com or from your local sales representative.

Application Areas

- Fabrication of temporary crowns, bridges, inlays, onlays, and veneers
- Fabrication of long-lasting temporary restorations
- Use of the material for prefabricated temporary crowns made of composite (e.g., Protemp™ Crown) and metal (e.g., Iso-Form crowns)

Precautionary Measures

For Patients

This product contains substances that may cause an allergic reaction in certain individuals. Avoid contact with this product in patients with known acrylate and/or peroxide allergies. If prolonged contact with soft tissue occurs, flush with large amounts of water. If allergic reactions occur, seek medical attention as needed, remove the product if necessary, and discontinue future use.

For Dental Personnel

This product contains substances that may cause an allergic reaction in certain individuals. Avoid contact with this product in patients with known acrylate and/or peroxide allergies. If prolonged contact with soft tissue occurs, flush with large amounts of water. If allergic reactions occur, seek medical attention as needed, remove the product if necessary, and discontinue future use.

Preparatory Work

- Remove the impression, silicone, or polyether impression from the intact row of teeth.
- Instead of an impression, a laboratory-made vacuum-formed template or, for single crowns, a preformed crown model may be used for shaping.
- To enhance the strength of the temporary restoration, remove interproximal grooves in the impression. If required due to large amounts of wax, enlarge the contact sites by cutting.
- For easy replacement in the mouth, relief undercuts.
- With missing teeth or gaps in the molars area, cut a groove into the impression to produce a stable, bar-shaped connection in the corresponding restoration.
- If required, prior to taking the impression, close any gaps between front teeth with denture teeth acting as spacers; interlock multiple denture teeth with wax.

Times

- 00:00–00:40 min:sec Place Paradigm™ Temporization Material into the impression or into the matrix/mold and reposition in the mouth
- 00:40–01:40 min:sec Setting in the mouth or in the matrix/mold (fringing)
- 01:40–02:50 min:sec Best point in time for removal of the impression or of the matrix/mold
- Up to 05:00 min:sec Complete setting
- From 05:00 min:sec Removal of the temporary restoration from the impression or from the matrix/mold (fringing)

The above processing times are applicable to products which have been stored and processed at a temperature of 23°C/72°F and at 50% relative humidity. Higher temperatures shorten the setting times and reduce the retention. The size of the restoration also has an impact on the setting times: as the size of the workpiece increases, the setting accelerates.

Dosing and Mixing

- The pastes are dosed and statically mixed in the Garant™ Dispenser.
 - Check whether both cartridge openings are free from obstruction.
 - Use the enclosed original Garant™ Mixing Tips Blue only!
- During the initial dosing**
- Remove and discard the cartridge lid. Under no circumstances should the cartridge lid be used to remove the cartridge as this could result in the inclusion of air bubbles.
 - Check whether both cartridge openings are free from obstruction and remove obstructing material, if any, with a suitable instrument.
 - Before a mixing tip is applied for the first time, squeeze paste out onto a mixing pad until the barrel and catalyst pastes are at the same level. Discard the paste which has been squeezed out.

Application

- Apply the material a certain amount of time to emerge from the newly attached mixing tip. The flow of material is interrupted as soon as the pressure on the handle is reduced. Material that has hardened during the mixing tip stay force, as this may damage both the cartridge and the mixing tip.
- Core build-ups which are fabricated using composites, composites, or resin-modified glass-ionomer cements must be sealed before the application of Paradigm™ Temporization Material, e.g., using Vaseline.
- First, install a new mixing tip.
- Prior to each application, extrude a small (pea-sized) amount of paste from the mixing tip onto a mixing pad and discard this material.

- Subsequently load the dried impression/matrix or vacuum-formed template in the relevant spaces from the bottom up.
- Reposition the impression or vacuum-formed template in the mouth/Mod on the model.
- The material attains a hard-elastic consistency within 1Min 40 sec after the onset of mixing. The material and the Impression or vacuum-formed template must be removed from the mouth/Mod within 2 min 30 sec after the onset of mixing.
- Check the setting process by observing the excess material in the mouth or on the matrix/mold.
- Keep the filled mixing tip on the cartridge to serve as a closure until the next use.

Finishing

- Completely remove the oxygen inhibition layer using alcohol (e.g., ethanol) once the material is completely cured (not earlier than 5 minutes after the onset of mixing). Polish if desired, e.g., with SoLo-lex™ polishing discs.
- After the product has cured completely, remove any residues.
- If necessary, finish and polish the temporary restoration using suitable carbide burs.
- Do not breathe in polishing dust. Use mouth protective device, safety glasses, and aspiration.

Individualization/Shape Correction

When using a Filtek Flowable Composite and Paradigm™ Temporization Material for custom-shaping and shape correction, the use of this material is no need for a bonding agent. If you intend to use other products, please refer to the pertinent manufacturer's Instructions for Use.

Compatible Products

- Flowable composites, e.g., from the Filtek™ Flowable Composite product line (e.g., A2, A3)
- Please comply with the corresponding Instructions for Use!
- Paradigm™ Temporization Material
- Processing: please refer to the information provided under "Application" and "Finishing"

Application of a Filtek Flowable Composite to New Temporary Restorations

In newly manufactured temporary restorations, the Filtek Flowable Composite can be applied either to the finished surface or to a surface which has been prepared by sanding or polishing in either case.

- Remove any contamination – e.g., from saliva or polishing dust – from the dispenser until a mix is ready.
- It is not necessary if the inhibition layer is clean or if the temporary restoration is finished and clean.
- For reasons of hygiene, dose the Filtek Flowable Composite into the dispenser until a mix is ready.
- Apply the composite to the temporary restoration in layers of maximum 1mm using a suitable instrument.
- Light-cure each layer separately for 20 sec.
- Trim the composite with carbide burs and polish, e.g., with SoLo-lex™ polishing discs.

Application of a Filtek Flowable Composite to Temporary Restorations That Were in Use

- Thoroughly roughen the surface of the temporary restoration with rotary instruments.
- Remove any contamination – e.g., from saliva or polishing dust – with alcohol, then dry the restoration with a stream of air.
- For reasons of hygiene, dose the Filtek Flowable Composite into the dispenser until a mix is ready.
- Apply the composite to the temporary restoration in layers of maximum 1mm using a suitable instrument.
- Light-cure each layer separately for 20 sec.
- Trim the composite with carbide burs and polish, e.g., with SoLo-lex™ polishing discs.

Lining of Prefabricated, Temporary Crowns Made of Composite and Metal

Core build-ups which are fabricated using composites, composites, or resin-modified glass-ionomer cements must be sealed before the application of Paradigm™ Temporization Material, e.g., using Vaseline.

Shape and cure prefabricated, temporary composite crowns (e.g., Protemp™ Crown) as follows:

- Please comply with the corresponding Instructions for Use!
- Remove any contamination – e.g., from saliva or polishing dust – with alcohol (e.g., ethanol), then dry with a stream of air.

Apply Paradigm™ Temporization Material to the prefabricated, temporary crown; see "Application".

- Place the prefabricated crown on the prepared tooth/die.
- The material attains a hard-elastic consistency within 1Min 40 sec after the onset of mixing and must be removed, along with the prefabricated crown, from the mouth/Mod within 2 min 30 sec after the onset of mixing.
- Check the setting process by observing the residue in the mouth or in the matrix/mold.
- Keep the filled mixing tip on the cartridge to serve as a closure until the next use.
- Remove the oxygen inhibition layer and finish the residues, see "Finishing".

Cementation

- If the prepared tooth or restoration has been sealed off (e.g., using Vaseline), remove the rest of the sealant thoroughly from the temporary restoration.
- Insert the temporary restoration with a commercial temporary cement, e.g., ReLayX™ Temp E or ReLayX™ Temp Neo. If the application is used as a long-lasting temporary restoration, we recommend using a strong adhesive such as ReLayX™ Unisem 2.
- Eugenol-containing cements may impair the setting of composite Luting cements to be used at a later time in the final luting step. If a composite cement is to be used for the final luting step, seal the temporary restoration with a eugenol-free cement.

After Cementation

- Thoroughly examine and explore the sulcus and surrounding dentition. Remove any residuals of the temporary material or cement from the mouth.
- Fabrication of a Temporary Restoration in the Laboratory Model Preparation**
- Take an impression to prepare an unsaved work model. If necessary, use a wax-up (wax modeling). Prepare a silicone matrix using a laboratory silicone.
 - Remove the matrix and wax from the model, e.g., using a steam cleaner or boiling water.
 - Use the use of a suitable dental wax to block out the undercuts on the work model will make the removal of the temporary restoration from the model easier. Apply a thin layer of wax to the die in addition. As necessary, cut out interproximal grooves and constrictions in the gingiva area in the matrix.
 - Trim the model with a suitable water-cooled coarse isolation in accordance with the applicable instructions for use.

Preparation of the Temporary Restoration

- Install a new mixing tip.
- Prior to each application, extrude a small (pea-sized) amount of paste from the mixing tip onto a mixing pad and discard this material.
- Subsequently load the dried impression/matrix in the relevant spaces from the bottom up.
- Reposition the impression/matrix or the vacuum-formed template in the mouth/Mod on the model.
- The material attains a hard-elastic consistency within 1Min 40 sec after the onset of mixing. The material and the Impression or vacuum-formed template must be removed from the mouth/Mod within 2 min 30 sec after the onset of mixing.
- Check the setting process by observing the excess material in the mouth or on the matrix/mold.
- Do not wait until the temporary restoration has hardened completely before removing it (see "Times").

- Keep the filled mixing tip on the cartridge to serve as a closure until the next use.

Finishing the Temporary Restoration

- To obtain an optimal composite surface, leave the temporary restoration in the matrix for at least 5 min (after start of mixing), then remove the matrix.
- Use alcohol (e.g., ethanol) to completely remove the oxygen inhibition layer in order to obtain a smooth, shiny surface.

Paradigm Temporization Material does not require any additional curing methods, they would cause excessive shrinkage. **Do not use heat, pressure, or light for any additional curing of the temporary restoration.**

Information about finishing and polishing can also be found in "Finishing" and "Individualization"/"Shape Correction". **Do not clean the temporary restoration with steam** because this can cause the material to become brittle. Do not use alcohol (e.g., ethanol) or disinfectants containing alcohol to clean the temporary restoration.

Note: Colored occlusal spray can cause discoloration of the temporary restoration.

Storage and Transport

- Ensure optimal fit of the temporary restoration by placing it on the model and transporting it on the model.

Repair

- Provide for mechanical and chemical retention by roughening the break surfaces and adjacent areas.
- Clean with alcohol and dry with air.
- Apply together with Paradigm™ Temporization Material or a Filtek Flowable Composite
- Cure Paradigm™ Temporization Material or the Filtek Flowable Composite as described in the relevant Instructions for Use.
- Once the material is fully set, finish as usual. When repaired with Paradigm™ Temporization Material, wait 10–15 min before loading.

Cleaning and Disinfection of the Garant™ cartridge

Warnings/Cautions/Dangers

- Make sure that the cleaning and disinfectant agents you have chosen do not contain any of the following materials:
 - Oxidation agents (e.g., hydrogen peroxides)
 - Oils
 - Glutaraldehydes
- Use only validated methods for cleaning and disinfection. Use solely cleaning and disinfectant agents tested for efficacy and compatibility with the medical device (e.g., FDA certification). Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals.

Limitations on Reprocessing

The medical device may be reused up to 500 times, provided that it is damaged and clean.

General

This medical device is non-sterile when delivered and must be cleaned and disinfected before initial use and prior to every further use.

Preparation for Decontamination

The medical device should be reprocessed as soon as possible after every use.

Reprocessing

The reprocessing procedure must be conducted before the initial use of the medical device before its initial use and after every further use of the medical device. The wipes contain the disinfectant Cavicide™ [alcohol-quaternary ammonium disinfectant]. Follow the cleaning agent manufacturer's directions for use.

Manual Cleaning by Wiping

- Ready-to-use disinfectant wipes (Cav Wipes™) are used to clean the medical device before its initial use and after every further use of the medical device. The wipes contain the disinfectant Cavicide™ [alcohol-quaternary ammonium disinfectant]. Follow the cleaning agent manufacturer's directions for use.
- Continue wiping the device with the wipes until no contamination is visible.
- Special care must be exercised in cleaning the grooves.

Manual Disinfection by Wiping (Intermediate Level)

- The device is disinfected with ready-to-use disinfectant wipes (Cav Wipes™). The wipes contain the disinfectant Cavicide™ [alcohol-quaternary ammonium disinfectant]. Follow the cleaning agent manufacturer's directions for use, observed in particular the contact times.
- Disinfect the medical device by using at least one new wipe for every wipe to wipe the device for a period of three minutes, ensuring that all of the surfaces of the medical device remain moist during the entire three-minute period.
- Special care must be exercised in disinfecting the grooves.

Inspection, Maintenance, Testing

- Examine the medical device for damage, discoloration, and deformation before every use.
- Do not use damaged devices under any circumstances.
- If the medical device is not visibly clean, repeat the reprocessing procedure.

Cleaning of the Garant Dispenser

- Remove the Garant™ cartridge.
- Remove paste that has not set using an alcohol-soaked tissue.
- Use the dispenser handle and plunger can be disassembled up to a temperature of 135°C/275°F; remove the plunger beforehand.
- Glutaraldehyde-based solutions can be used for disinfection.

Storage and Stability

Store the product at 15–25°C/59–77°F. Do not store in the refrigerator. Do not use after the expiration date.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH does not warrant THE WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for use. If this material is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the product as described by 3M Deutschland GmbH product.

Limitation of Liability

Notwithstanding to the extent provided by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether for direct, indirect, special, incidental or consequential, regard of the tortious asserted, including warranty, contract, negligence or strict liability.

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3M, ESPE, Filtek, Garant, Penta, Protemp, ReLayX and SoLo-lex are trademarks of 3M or 3M Deutschland GmbH.
© 2016, 3M. All rights reserved.